



$$Q \propto \sqrt{\frac{(P_{G1}-P_{G2})}{(P_{G1}+P_{G2})+2*Patmosphere}}$$

An important and necessary simplification is that P_{G1} and P_{G2} encountered in our system are of the order of tens of cmH₂O while $Patmosphere$ is of the order of a thousand cmH₂O of pressure. At sea level, $Patmosphere$ is approximately 1000 cmH₂O. Even at an altitude of 15,000 feet, $Patmosphere$ is approximately 600 cmH₂O. On the other hand, the P_{G1} and P_{G2} range from 1 cmH₂O to 60 cmH₂O.

Thus the $(P_{G1} + P_{G2})$ term is negligible compared to $(2 * Patmosphere)$ even more so since it is preceded by a square root. The flow equation can be simplified to the one below.

$$Q \propto \sqrt{\frac{(P_{G1}-P_{G2})}{Patmosphere}}$$

Recalling *Equation 2* from the theory above, this equation can be recast as below given that the orifice characteristics and pressure variations are the same for every system.

$$Q = C * \frac{\sqrt{(P_{G1}-P_{G2})}}{\sqrt{(Patmosphere)}} \text{ where } C = f(Re) \text{ Reynold's number}$$

At a given point in time, $Patmosphere$ is also a constant. So, the above equation further reduces to the one below.

$$Q = \left(\frac{C}{\sqrt{(Patmosphere)}} \right) * \sqrt{(P_{G1}-P_{G2})}$$

The equation needs further simplification to ease the computation burden of the square root computation for and inexpensive micro-controller. The constraints are as below.

A Simplified Emergency Respiration Assist Device



INSPIRE-100

Salient Features

*Emergency Adult
Respiratory
Support*

*Support from
Initiation to
Weaning*

*No Need for
Compressed air*

*Connect to O₂
Cylinder or
Concentrator*

*Full range of
Respiration
Parameters*

*Breath
Synchronization
Patient Comfort*

*Elaborate
Remote
Monitoring*

*Single Limb
Breathing Circuit*

*Compatible with
Standard
Accessories*

*Easy-to-use
Easy-to-train
Rugged & Robust*

*Power
Consumption
100W*

*Eminently
Budget Friendly*

Breathing System Components



Commonly used Ventilation Modes

CMV	Continuous Mandatory Ventilation
ACV	Synchronized Assist Control Ventilation
SIMV	Synchronized Intermittent Mandatory Ventilation
PSV	Pressure Support Ventilation

Full Set of Alarm Alerts

Max Pressure	Pressure Leak	Pressure Loss
Airway Blockage	System Temperature	Coughing Hiccupping
Inconsistent Parameters	Extreme Parameter Combination	And many more ...





Volume Controlled Breaths

Tidal Volume 200 to 600 ml	Respiratory Rate 10 to 30 bpm	I:E Ratio 1:1 1:2 1:3
PEEP 4 to 15 cmH ₂ O	Max Pressure 15 to 60 cmH ₂ O	FiO₂ Support System Managed Externally Controlled


Pressure Supported Breaths

Support Pressure	5 cmH ₂ O to 35 cmH ₂ O
Support Pressure Termination	Flow Triggered Time Triggered


Elaborate Remote Monitoring




Menu-driven, Tactile, Intuitive,
Easy-to-read Control Panel



Live Dashboard
Recording, Analyzer



Live Dashboard for
Nurses' Station



Field Upgradeable

